


Tc Reconstruction

User Manual



(See Fig.S2, main text; yellow: SMLM&tcSMLM, red: SMLM; green: tcSMLM)

Table of Contents

Getting Started	3
1. Requirements:	3
2. Installation:	3
3. Main scripts and examples	3
Example (1):	4
Example (2):	4
Example (3):	4
4. Finding the optimal MW	4
5. Final Reconstruction	6

Getting Started

Here, we provide a set of instructions toward employing tcReconstruction.

1. Requirements:

MATLAB. Has been tested on R2016b but scripted to work on all versions from 2014-onward.

Parallel Computing Tollbox:

<https://www.mathworks.com/products/parallel-computing.html>

2. Installation:

A. Download the folder 'tcReconstruction' from our GitHub account:

<https://github.com/ShermanLab/tcSMLM> and unzip it.

B. Choose 'tcReconstruction' as the main folder in your MATLAB.

C. Open "examples.mat" in MATLAB.

D. Download the following samples:

https://drive.google.com/file/d/1p_OeiV6NGrte2murSd4b5LyteKR-sTrL/view?usp=sharing

<https://drive.google.com/file/d/1eYBUwUtZwMwMh1c0-VyRX1DkHIVNbjiJ/view?usp=sharing>

Insert the samples into the tcReconstruction folder.

Now, in the workspace, there are two RawData sets:

- 1.) Alexa647 experimental data with one emitter (See methods)
- 2.) Simulation data with two emitters.

In addition, there is one .nd2 file in your current folder, one .tif file and the tcReconstruction functions.

3. Main scripts and examples

- A. tcData is created using the matlab script named "tcData.m". In order to apply the script, a 3D matlab matrix containing the RawData should be available (either

created using dSTORM/STORM method or by other methods containing temporal information). The tcData function has 3 parameters:

time_step: time_step = 1 is recommended;

MW: see below how to choose the optimal MW;

RawData: the original, unresolved data (intensity information).

Example (1):

```
tc_Alexa647 = tcData(1,65,Alexa647);
```

- B. In case you use an .nd2 file, a matlab function named “Import_nd2.m” is available. The function also creates a .tif stack file in the folder containing the .nd2 file.

Example (2):

```
PAmC_NRas = Import_nd2('Desktop/tcReconstruction', 'PAmC_NRas', 2000, 0);
```

- C. Otherwise, if you have a RawData in a .tif stack file, a simple function named “ReadTif.m” is available. The function provides a faster way to import .tif files comparable to a direct import using Matlab.

Example (3):

```
PAmC_NRas = ReadTif('PAmC_NRas.tif', 2000);
```

4. Finding the optimal MW

The default parameters are in the range: $25 \leq MW \leq 55$. Lower for slow frame rate (20-50 fps), higher for fast frame rate (100+ fps). The function to find the optimal MW is: ‘FindBestMovingWindow.m’. The function has 4 parameters:


X = an array of size n.

Y = an array of size n.

Smooth = moving average. Ranges between 10-20, depends on the noise in the data.

rawdata = the original (experimental, as intensity) data.

For example, we look at the “PAmC_NRas” sample. Here is the sum intensity where we chose a number of locations:



These locations contain the X and Y arrays:

$X = [40, 112, 14, 28, 21, 57, 86]$, $Y = [16, 29, 112, 124, 137, 132, 107]$.

So:

`FindBestMovingWindow([40, 112, 14, 28, 21, 57, 86], [16, 29, 112, 124, 137, 132, 107], 20, PAmC_NRas);`

We receive 3 figs in return. The important of them is the third one, containing the MW for each curve. Here, we find: the minimal MW which is bigger that 20 (see main text, Figs5&6) is $MW = 27$. Therefore, the optimal tcData for PAmC_NRas is:

```
tc_PAmC_NRas = tcData(1,27, PAmC_NRas);
```

5. Final Reconstruction

Once the *tcData* is created any SMLM / super-resolution algorithm could be applied. That could be done by saving the data to an external data file (such as .tif) or using a MATLAB based super resolution algorithm.